

NATIONAL STANDARD METHOD

INVESTIGATION OF GASTRIC BIOPSIES FOR *HELICOBACTER PYLORI*

BSOP 55

Issued by Standards Unit, Evaluations and Standards Laboratory
Specialist and Reference Microbiology Division

Association of Medical Microbiologists
Association of Medical Microbiologists



INVESTIGATION OF GASTRIC BIOPSIES FOR *HELICOBACTER PYLORI*

Issue no: 4.1 Issue date: 03.05.05 Issued by: Standards Unit, Evaluations and Standards Laboratory Page 1 of 14
Reference No: BSOP 55i4.1

This SOP should be used in conjunction with the series of other SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

STATUS OF NATIONAL STANDARD METHODS

National Standard Methods, which include standard operating procedures (SOPs), algorithms and guidance notes, promote high quality practices and help to assure the comparability of diagnostic information obtained in different laboratories. This in turn facilitates standardisation of surveillance underpinned by research, development and audit and promotes public health and patient confidence in their healthcare services. The methods are well referenced and represent a good minimum standard for clinical and public health microbiology. However, in using National Standard Methods, laboratories should take account of local requirements and may need to undertake additional investigations. The methods also provide a reference point for method development.

National Standard Methods are developed, reviewed and updated through an open and wide consultation process where the views of all participants are considered and the resulting documents reflect the majority agreement of contributors.

Representatives of several professional organisations, including those whose logos appear on the front cover, are members of the working groups which develop National Standard Methods. Inclusion of an organisation's logo on the front cover implies support for the objectives and process of preparing standard methods. The representatives participate in the development of the National Standard Methods but their views are not necessarily those of the entire organisation of which they are a member. The current list of participating organisations can be obtained by emailing standards@hpa.org.uk.

The performance of standard methods depends on the quality of reagents, equipment, commercial and in-house test procedures. Laboratories should ensure that these have been validated and shown to be fit for purpose. Internal and external quality assurance procedures should also be in place.

Whereas every care has been taken in the preparation of this publication, the Health Protection Agency or any supporting organisation cannot be responsible for the accuracy of any statement or representation made or the consequences arising from the use of or alteration to any information contained in it. These procedures are intended solely as a general resource for practising professionals in the field, operating in the UK, and specialist advice should be obtained where necessary. If you make any changes to this publication, it must be made clear where changes have been made to the original document. The Health Protection Agency (HPA) should at all times be acknowledged.

The HPA is an independent organisation dedicated to protecting people's health. It brings together the expertise formerly in a number of official organisations. More information about the HPA can be found at www.hpa.org.uk.

The HPA aims to be a fully Caldicott compliant organisation. It seeks to take every possible precaution to prevent unauthorised disclosure of patient details and to ensure that patient-related records are kept under secure conditions¹.

More details can be found on the website at www.evaluations-standards.org.uk. Contributions to the development of the documents can be made by contacting standards@hpa.org.uk.

Please note the references are now formatted using Reference Manager software. If you alter or delete text without Reference Manager installed on your computer, the references will not be updated automatically.

Suggested citation for this document:

Health Protection Agency (2004). *Investigation of gastric biopsies for Helicobacter pylori*. National Standard Method BSOP 55 Issue 4. http://www.hpa-standardmethods.org.uk/pdf_bacteriology.asp.

INVESTIGATION OF GASTRIC BIOPSIES FOR *HELICOBACTER PYLORI*

Issue no: 4.1 Issue date: 03.05.05 Issued by: Standards Unit, Evaluations and Standards Laboratory Page 2 of 14

Reference No: BSOP 55i4.1

This SOP should be used in conjunction with the series of other SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

INDEX

- STATUS OF NATIONAL STANDARD METHODS..... 2**
- INDEX..... 3**
- AMENDMENT PROCEDURE 4**
- SCOPE OF DOCUMENT 5**
- INTRODUCTION..... 5**
 - TECHNICAL INFORMATION 6
- 1.0 SAFETY CONSIDERATIONS..... 8**
 - 1.1 SPECIMEN COLLECTION 8
 - 1.2 SPECIMEN TRANSPORT AND STORAGE 8
 - 1.3 SPECIMEN PROCESSING..... 8
- 2.0 SPECIMEN COLLECTION 8**
 - 2.1 OPTIMAL TIME OF SPECIMEN COLLECTION 8
 - 2.2 CORRECT SPECIMEN TYPE AND METHOD OF COLLECTION 8
 - 2.3 ADEQUATE QUANTITY AND APPROPRIATE NUMBER OF SPECIMENS 8
- 3.0 SPECIMEN TRANSPORT AND STORAGE..... 8**
 - 3.1 TIME BETWEEN SPECIMEN COLLECTION AND PROCESSING..... 8
 - 3.2 SPECIAL CONSIDERATIONS TO MINIMISE DETERIORATION 8
- 4.0 SPECIMEN PROCESSING..... 9**
 - 4.1 TEST SELECTION 9
 - 4.2 APPEARANCE 9
 - 4.3 MICROSCOPY 9
 - 4.4 CULTURE AND INVESTIGATION 9
 - 4.5 IDENTIFICATION 10
 - 4.6 ANTIBIOTIC SUSCEPTIBILITY TESTING 10
- 5.0 REPORTING PROCEDURE 10**
 - 5.1 MICROSCOPY 10
 - 5.2 CULTURE 11
 - 5.2 ANTIMICROBIAL SUSCEPTIBILITY TESTING 11
- 6.0 REPORTING TO THE HPA (LOCAL AND REGIONAL SERVICES AND CDSC CENTRE FOR INFECTIONS)..... 11**
- REFERENCES 13**

INVESTIGATION OF GASTRIC BIOPSIES FOR *HELICOBACTER PYLORI*

Issue no: 4.1 Issue date: 03.05.05 Issued by: Standards Unit, Evaluations and Standards Laboratory Page 3 of 14

Reference No: BSOP 55i4.1

This SOP should be used in conjunction with the series of other SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

AMENDMENT PROCEDURE

Controlled document reference	BSOP 55
Controlled document title	Standard Operating Procedure for the investigation of gastric biopsies for <i>Helicobacter pylori</i>

Each National Standard Method has an individual record of amendments. The current amendments are listed on this page. The amendment history is available from standards@hpa.org.uk.

On issue of revised or new pages each controlled document should be updated by the copyholder in the laboratory.

Amendment Number/ Date	Issue no. Discarded	Insert Issue no.	Page	Section(s) involved	Amendment
5/ 03.05.05	4	4.1	1	Front page	Redesigned
			2	Status of document	Reworded
			4	Amendment page	Redesigned

INVESTIGATION OF GASTRIC BIOPSIES FOR *HELICOBACTER PYLORI*

Issue no: 4.1 Issue date: 03.05.05 Issued by: Standards Unit, Evaluations and Standards Laboratory Page 4 of 14
Reference No: BSOP 55i4.1

This SOP should be used in conjunction with the series of other SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

STANDARD OPERATING PROCEDURE FOR THE INVESTIGATION OF GASTRIC BIOPSIES FOR *HELICOBACTER PYLORI*

Type of specimen: Gastric biopsy

SCOPE OF DOCUMENT

This SOP describes the processing and bacteriological investigation of gastric biopsies for *Helicobacter pylori*.

INTRODUCTION

Infection with *H. pylori* is associated with peptic ulceration and is a risk factor for gastric cancer. There is increasing evidence that it plays an important role in non-ulcer dyspepsia. Acute symptoms of gastritis and epigastric pain, nausea and vomiting may occur and usually subside, but hyperchlorhydria may persist for much longer².

Gastric biopsies are the specimens of choice for the culture of *H. pylori*, but culture methods have been applied to other specimens³.

Diagnosis of *H. pylori* infection may be established by a number of invasive and non-invasive techniques⁴.

Invasive techniques examine gastric biopsies taken at endoscopy and include⁴⁻⁶:

- culture of the organism
- histology
- biopsy urease test
- microscopy
- polymerase chain reaction (PCR)

Culture of the organism is the most specific method and offers opportunity for conventional antimicrobial susceptibility testing if required. This is important in predicting and evaluating the efficacy of treatment, and in identifying re-infections.

Histological examination is as sensitive as culture for the detection of *H. pylori*, and has a high degree of specificity⁷.

Neither culture nor histology will provide a rapid diagnosis.

The biopsy urease test is rapid, sensitive and cost effective. Positive results are often available within minutes but negative reporting may take a great deal longer, according to manufacturers instructions. It is recommended for use in combination with either culture or histology, depending on local facilities. This test is often carried out in the endoscopy suite. Commercial kits are available which are highly accurate but are also expensive.

INVESTIGATION OF GASTRIC BIOPSIES FOR *HELICOBACTER PYLORI*

Issue no: 4.1 Issue date: 03.05.05 Issued by: Standards Unit, Evaluations and Standards Laboratory Page 5 of 14

Reference No: BSOP 55i4.1

This SOP should be used in conjunction with the series of other SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

Microscopy of tissue smears - organisms may be stained using Giemsa or Gram stains according to preference. Sensitivities of up to 90% have been reported if two biopsies are examined, but this method requires technical expertise. It is the only other rapid method than the biopsy urease test.

PCR has been used for the detection of *H. pylori* in various samples, although its role in routine diagnosis remains to be established.

Non-invasive techniques (avoiding the need for expensive and invasive endoscopy) include:

- serology
- urea breath tests (UBTs)
- faecal antigen tests

Serology ELISA based tests have become more accurate in recent years. A large number of commercial kits are available. IgG detection is commonly used and has the greatest published evidence, but IgA is also available and inflammation marker are sometimes included. Several point-of-care tests using whole blood are available.

False positives have been shown to increase with the age of the patient⁸. The use of this technique to confirm eradication is limited by the variable but prolonged presence of immunoglobulins post clearance.

Urea breath tests (UBT's) are considered to be the diagnostic gold standard⁹. The urea molecule used is labelled with either ¹⁴C or ¹³C. The former employs simple instrumentation but the radioactive nature of the test inhibits its use. The latter uses a stable isotope but requires complex instrumentation. Several ¹³C labelled tests are available commercially as postal services or using dedicated in-house instrumentation. Local methods can be created if the laboratory has access to mass spectrometry.

UBT's allow the rapid assessment of eradication efficacy however sample collection requires time and technical understanding.

Faecal Antigen Tests are the newest development in the range of diagnostic options. Little independent information has been published but what has been indicates that the test is sensitive and specific¹⁰ and has the ability to confirm eradication.

These tests can be performed remotely from the patient allowing batched processing.

TECHNICAL INFORMATION

Optimal growth requirements for the isolation of *H. pylori* are a moist, microaerobic atmosphere of 5-7% O₂ and 5-10% CO₂ at 35-37°C^{3,4}. Gas generating kits for microaerobic conditions are commercially available.

Homogenisation may be performed, but it is more time consuming and requires the use of a Griffiths grinder or an unbreakable alternative. Biopsy can be cut finely using sterile scalpel.

Cultures should be incubated for up to 7 days, although colonies are usually visible at 3-5 days⁴.

Gram stain using dilute carbol fuchsin counterstain enhances morphology.

There is currently no single medium that is best for the isolation of *H. pylori* although blood based media are preferred. Several have been described^{4,11-13}.

Blood-free media, containing alternative supplements, may not be as good for primary isolation¹².

INVESTIGATION OF GASTRIC BIOPSIES FOR *HELICOBACTER PYLORI*

Issue no: 4.1 Issue date: 03.05.05 Issued by: Standards Unit, Evaluations and Standards Laboratory Page 6 of 14

Reference No: BSOP 55i4.1

This SOP should be used in conjunction with the series of other SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

Antimicrobial supplements may be added to media to inhibit overgrowth with contaminating bacteria and fungi¹⁴.

H. pylori is sensitive to clindamycin, cephalosporins and sodium desoxycholate, none of which should be used in the selective medium¹². Selective media for *Neisseria gonorrhoeae* may be used although about 5% of isolates of *H. pylori* may be inhibited by colistin or polymixin B contained in the medium.

Contamination with moulds may be reduced by the incorporation of an antifungal agent into the media such as cyclohexamide (100mg/L)³ and thorough cleaning of equipment before and after use. Autoclaving of jars previously contaminated with moulds is recommended.

Best results are obtained using both selective and non-selective media¹⁵.

Confirmation of the organism relies on the characteristic "seagull" morphology in the Gram film, and positive oxidase and rapid urease tests.

INVESTIGATION OF GASTRIC BIOPSIES FOR *HELICOBACTER PYLORI*

Issue no: 4.1 Issue date: 03.05.05 Issued by: Standards Unit, Evaluations and Standards Laboratory Page 7 of 14

Reference No: BSOP 55i4.1

This SOP should be used in conjunction with the series of other SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

1.0 SAFETY CONSIDERATIONS¹⁶⁻²⁶

1.1 SPECIMEN COLLECTION

N/A

1.2 SPECIMEN TRANSPORT AND STORAGE

Sterile leakproof container in a sealed plastic bag

1.3 SPECIMEN PROCESSING

Containment Level 2

The above guidance should be supplemented with local COSHH and risk assessments

2.0 SPECIMEN COLLECTION

2.1 OPTIMAL TIME OF SPECIMEN COLLECTION

Before antimicrobial therapy where possible

2.2 CORRECT SPECIMEN TYPE AND METHOD OF COLLECTION

Gastric biopsy specimens are usually taken from the gastric antrum at endoscopy, sometimes from the body depending on location of inflammation

2.3 ADEQUATE QUANTITY AND APPROPRIATE NUMBER OF SPECIMENS

At the discretion of the endoscopist as it depends on the individual patient

3.0 SPECIMEN TRANSPORT AND STORAGE

3.1 TIME BETWEEN SPECIMEN COLLECTION AND PROCESSING

Specimens should be transported and processed as soon as possible (preferably within 6h)⁴

3.2 SPECIAL CONSIDERATIONS TO MINIMISE DETERIORATION

It is important to maintain a moist atmosphere during transport

Where culture is to be carried out within 6h⁴:

The biopsy should be placed in a small, sterile container such as a bijou bottle, containing a small amount (approximately 100µL) of sterile isotonic saline to preserve moisture. Dent's transport medium can be used¹⁴. This is available from the Helicobacter Reference Unit

NOTE: Sensitivity of the microscopy may be reduced if the biopsy is submerged in the saline, mucus globules form and production of a satisfactory smear becomes difficult

INVESTIGATION OF GASTRIC BIOPSIES FOR *HELICOBACTER PYLORI*

Issue no: 4.1 Issue date: 03.05.05 Issued by: Standards Unit, Evaluations and Standards Laboratory Page 8 of 14

Reference No: BSOP 55i4.1

This SOP should be used in conjunction with the series of other SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

Where delays of >6h are expected^{4,15}:

The biopsy should be covered with approximately 1ml brain heart infusion broth in a small sterile container, such as a bijou bottle, and stored at 4°C for up to 48h Dent's transport medium can be used. Organisms will remain viable in Amies transport medium, but if this is used, care is required to ensure that the mucosa has not become detached from the rest of the biopsy¹⁴

Biopsies may be stored for up to 6 months at -70°C in broth containing 20-25% glycerol although viability will be reduced

4.0 SPECIMEN PROCESSING

4.1 TEST SELECTION

The biopsy urease test is often performed in the endoscopy suite so only culture and microscopy may be required in the laboratory

The order in which any/all of the tests are performed will be in accordance with local protocol

4.2 APPEARANCE

N/A

4.3 MICROSCOPY

4.3.1. Standard

Pick up the biopsy with a sterile swab and smear vigorously on to a clean microscope slide (a **sterile** slide is required if microscopy is performed before culture)

Staining and examination of the stained preparation need only be performed if the culture result is negative and the biopsy urease test positive. Gram or Giemsa stains are suitable

NOTE: Methods for staining procedures are contained in separate SOPs

4.4 CULTURE AND INVESTIGATION

4.4.1. Pre-treatment

N/A

4.4.2. Specimen processing

Culture

Using the same swab containing the biopsy as for microscopy (if performed), inoculate each agar plate (see BSOP 54)

For the isolation of individual colonies, spread inoculum using a sterile loop

NOTE: The simultaneous subculture of known control strains of *H. pylori* is recommended, especially if susceptibility testing is to be performed (see BSOP 45)

Metronidazole sensitive strain - NCTC 12822

INVESTIGATION OF GASTRIC BIOPSIES FOR *HELICOBACTER PYLORI*

Issue no: 4.1 Issue date: 03.05.05 Issued by: Standards Unit, Evaluations and Standards Laboratory Page 9 of 14

Reference No: BSOP 55i4.1

This SOP should be used in conjunction with the series of other SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

Metronidazole resistant strain - NCTC 12823

Biopsy urease test

Squash the biopsy on the end of the swab into urease broth

The swab should be broken off into the broth and left *in situ* throughout the test

Incubate the urease broth at ambient temperature for up to 24h

4.4.3 Culture media, conditions and organisms

For all specimens:

Clinical details/ conditions	Standard media	Incubation			Cultures read	Target organism
		Temp °C	Atmos	Time		
Gastritis Gastric biopsy	<i>H. pylori</i> selective agar	35-37	microaerobic moist chamber	7d	at 4-5d and 7d	<i>H. pylori</i>
	Blood agar	35-37	microaerobic moist chamber	7d	at 4-5d and 7d	

For these situations, add the following:

Clinical details/ conditions	Supplementary media	Incubation			Cultures read	Target organism
		Temp °C	Atmos	Time		
Biopsy urease test if not already performed in endoscopy suite	Biopsy urease broth	ambient	air	24h	hourly up to 6h and again at 24h	<i>H. pylori</i>

4.5 IDENTIFICATION

4.5.1. Minimum level in the laboratory

H. pylori: to species level

4.5.2. Referral to Reference Laboratories

Not routinely referred

Isolates with unusual or unexpected resistance, and whenever there is a laboratory or clinical problem or anomaly that requires elucidation

4.6 ANTIBIOTIC SUSCEPTIBILITY TESTING

Refer to SOP on Susceptibility Testing (BSOP 45)

5.0 REPORTING PROCEDURE

5.1 MICROSCOPY

Gram stain (if performed)

INVESTIGATION OF GASTRIC BIOPSIES FOR *HELICOBACTER PYLORI*

Issue no: 4.1 Issue date: 03.05.05 Issued by: Standards Unit, Evaluations and Standards Laboratory Page 10 of 14

Reference No: BSOP 55i4.1

This SOP should be used in conjunction with the series of other SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

Report presence or absence of *H. pylori*-like organisms

5.1.1. Microscopy reporting time

N/A

5.2 CULTURE

The following as appropriate:

Positive results: "*H. pylori* isolated"

Negative results: "*H. pylori* NOT isolated"

Biopsy urease test: Report biopsy urease test result as positive or negative

5.1.2. Culture reporting time

Clinically urgent culture results to be telephoned or sent electronically

Written report: 24h for biopsy urease test (if not already performed in the endoscopy suite), stating that a further report on the culture will be issued

Culture result within 7 days

Supplementary investigations: up to 7 days for microscopy

5.2 ANTIMICROBIAL SUSCEPTIBILITY TESTING

Report susceptibilities as clinically indicated

6.0 REPORTING TO THE HPA²⁷ (LOCAL AND REGIONAL SERVICES AND CDSC CENTRE FOR INFECTIONS)

Refer to the following:

Individual SOPs on organism identification

Health Protection Agency publications:

"Reporting to the CDR: A guide for laboratories"

"Hospital infection control: Guidance on the control of infection in hospitals"

Refer to current guidelines on CDSC and COSURV reporting

Local guidelines

INVESTIGATION OF GASTRIC BIOPSIES FOR *HELICOBACTER PYLORI*

Issue no: 4.1 Issue date: 03.05.05 Issued by: Standards Unit, Evaluations and Standards Laboratory Page 11 of 14


Reference No: BSOP 55i4.1

This SOP should be used in conjunction with the series of other SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

Appendix 1. TESTING ALGORITHM: GASTRIC BIOPSIES FOR *HELICOBACTER PYLORI*



INVESTIGATION OF GASTRIC BIOPSIES FOR *HELICOBACTER PYLORI*

Issue no: 4.1 Issue date: 03.05.05 Issued by: Standards Unit, Evaluations and Standards Laboratory Page 12 of 14

Reference No: BSOP 55i4.1

This SOP should be used in conjunction with the series of other SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

REFERENCES

1. Department of Health NHS Executive. The Caldicott Committee. Report on the review of patient-identifiable information. 1997. London.
2. Parsonnet J. *Helicobacter*. In: Gorbach SL, Bartlett JG, Blacklow NR, editors. *Infectious Diseases*. Philadelphia: WB Saunders Company; 1998. p. 1952-61.
3. Tompkins D. Microbiological tests. In: Northfield TC, Mendall M, Goggin PM, editors. *Helicobacter pylori* infection. Dordrecht: Kluwer Academic Publishers; 1993. p. 116-26.
4. Glupczynski Y. The diagnosis of *Helicobacter pylori* infection: a microbiologist's perspective. *Rev Med Microbiol* 1994;5:199-208.
5. Thillainayagam AV, Arvind AS, Cook RS, Harrison IG, Tabaqchali S, Farthing MJ. Diagnostic efficiency of an ultrarapid endoscopy room test for *Helicobacter pylori*. *GUT* 1991;32:467-9.
6. Goodwin CS, Armstrong JA. Microbiological aspects of *Helicobacter pylori* (*Campylobacter pylori*). *Eur J Clin Microbiol Infect Dis* 1990;9:1-13.
7. McNulty CA, Wise R. Rapid diagnosis of *Campylobacter*-associated gastritis. *Lancet* 1985;1:1443-4.
8. Schembri MA, Lin SK, Lambert JR. Comparison of commercial diagnostic tests for *Helicobacter pylori* antibodies. *J Clin Microbiol* 1993;31:2621-4.
9. Logan RP. Urea breath tests in the management of *Helicobacter pylori* infection. *Gut* 1998;43 Suppl 1:S47-S50.
10. Vaira D, Vakil N, Menegatti M, van't Hoff B, Ricci C, Gatta L, et al. The stool antigen test for detection of *Helicobacter pylori* after eradication therapy. *Ann Intern Med* 2002;136:280-7.
11. Hachem CY, Clarridge JE, Evans DG, Graham DY. Comparison of agar based media for primary isolation of *Helicobacter pylori*. *J Clin Pathol* 1995;48:714-6.
12. Tompkins DS. Isolation and characteristics of *Helicobacter pylori*. In: Rathbone B, Heatley RV, editors. *Helicobacter pylori* and infectious disease. Oxford: Blackwell Scientific Publications; 1992. p. 19-28.
13. Henriksen TH, Brorson O, Schoyen R, Thoresen T, Setegn D, Madebo T. Rapid growth of *Helicobacter pylori*. *Eur J Clin Microbiol Infect Dis* 1995;14:1008-11.
14. Dent JC, McNulty CA. Evaluation of a new selective medium for *Campylobacter pylori*. *Eur J Clin Microbiol Infect Dis* 1988;7:555-8.
15. Tompkins D. Diagnosis of *Helicobacter pylori* infection. *PHLS Microbiol Dig* 1997;14:34-6.
16. Advisory Committee on Dangerous Pathogens. Categorisation of biological agents according to hazard and categories of containment. 4th ed. Suffolk: HSE Books; 1995. Supplements 1, 1998 and 2, 2000.
17. Public Health Laboratory Service Standing Advisory Committee on Laboratory Safety. Safety Precautions: Notes for Guidance. 4th ed. London: Public Health Laboratory Service (PHLS); 1993.
18. Control of Substances Hazardous to Health Regulations 2002. General COSHH Approved Code of Practice and Guidance, L5. Suffolk: HSE Books; 2002.

INVESTIGATION OF GASTRIC BIOPSIES FOR *HELICOBACTER PYLORI*

Issue no: 4.1 Issue date: 03.05.05 Issued by: Standards Unit, Evaluations and Standards Laboratory Page 13 of 14

Reference No: BSOP 55i4.1

This SOP should be used in conjunction with the series of other SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

19. Health and Safety Executive. 5 steps to risk assessment: a step by step guide to a safer and healthier workplace, IND (G) 163 (REVL). Suffolk: HSE Books; 2002.
20. Health and Safety Executive. A guide to risk assessment requirements: common provisions in health and safety law, IND (G) 218 (L). Suffolk: HSE Books; 2002.
21. Health Services Advisory Committee. Safety in Health Service laboratories. Safe working and the prevention of infection in clinical laboratories and similar facilities. 2nd ed. Suffolk: HSE Books; 2003.
22. NHS Estates. Health Building Note 15. Accommodation for pathology services. 1st ed. London: Her Majesty's Stationary Office (HMSO); 1991. (Out of print - 2nd edition in press).
23. BS EN 12469: 2000. Biotechnology - performance criteria for microbiological safety cabinets. London: British Standards Institution (BSI); 2000.
24. BS 5726: 1992. Microbiological safety cabinets. Part 2. Recommendations for information to be exchanged between purchaser, vendor and installer and recommendations for installation. London: British Standards Institution (BSI); 1992.
25. BS 5726: 1992. Microbiological safety cabinets. Part 4. Recommendations for selection, use and maintenance. London: British Standards Institution (BSI); 1992.
26. Advisory Committee on Dangerous Pathogens. The management, design and operation of microbiological containment laboratories. Suffolk: HSE Books; 2001.
27. PHLS, CDSC. Reporting to the PHLS Communicable Disease Surveillance Centre: a reference for laboratories. May. 2001.

INVESTIGATION OF GASTRIC BIOPSIES FOR *HELICOBACTER PYLORI*

Issue no: 4.1 Issue date: 03.05.05 Issued by: Standards Unit, Evaluations and Standards Laboratory Page 14 of 14

Reference No: BSOP 55i4.1

This SOP should be used in conjunction with the series of other SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk